

8-Channel Variable Repetitive Syringe Dispenser

MODEL 8800

- Syringe type 8 channel dispenser.
- Suitable for dispensing to 96 well micro plates.
- Up to eight syringes can be selected for dispensing into various types of tissue culture plates.
- Six different volumes between 10 μ L – 100 μ L.
- The dispensing speed can be altered.

**Specifications**

Cat. No	Description	Q'ty
00-8800	MODEL 8800 (cartridge A&B one each, syringe A&B eight each)	1

Accuracy and Precision

syringe \ dial No	1	2	3	Accuracy(%)	Precision(%)
Syringe A (520μL)	10μL	20μL	40μL	±2.0	≤1.2
Syringe B (1300μL)	25μL	50μL	100μL		
Maximum number of dispenses	50	25	12	Size:285×82×23mm (with syringes set)	

Accessory

Accuracy _____

Stand _____

8-Channel Variable Repetitive Syringe Dispenser
MODEL 8800

Cat. No	Description	Q'ty
00-88-SAR	Syringe A (520 μ L)	100
00-88-SBR	Syringe B (1300 μ L)	100
00-88-CA	Cartridge A (for syringe A)	10
00-88-CB	Cartridge B (for syringe B)	10
00-88-TY	Tray	50
00-8800S	Storing capacity (Acrylic resin) (capacity:2pcs)	1

**Ergonomically Designed Repetitive Dispenser Nichimate Stepper**

Cat. No00-NS-STD
(Acrylic resin)
Storing capacity ...2pcs

**Variable Repetitive Syringe Dispenser MODEL 8100**

Cat. No00-8100S
(Acrylic resin)
Storing capacity ...2pcs